

New Products / Nuevo Productos



1575
Bar Holder (inox, Silicon) autoclave 160°
Fresero (fresas, contra angulo, pieza de mano)



1151
Cotton Holder (inox)



1154
Spirit Lamp (inox)





1145
Endodontic Box (inox)




1666
5 Pcs

Sterilize Trays INOX

1667
10 Pcs



Sterilize Trays Large INOX
1576



1441
Table of save clamps (inox)
porta para guardar clamps.
117x61x10 mm



1430
Led Curing Light

- Feature:
1. High articulation optical fiber (08 mm)
 2. Three working modes: strong, gradually strong, flashing
 3. Optional working time: 10s, 20s, 30s, 40s
 4. Light intensity: 1200mw/cm2
 5. Automatic inactive, automatic power off.



1429
Mixing Bowl Inox
30 mm